Oxygen CiTiceL® Specification



AO2 CiTiceL®

with Molex connector

Performance Characteristics

Output | 9 - 13mV in Air

Range 0-100% O₂

Resolution 0.01% O₂

Expected Operating Life | 360000%O,hrs at 20°C

286000%O₂hrs at 40°C

or 2 years in air at STP

T₉₀ Response Time | <5 seconds T₉₀₅ Response Time* | <40 seconds

Signal in 100%O₂ 100±1%

Linearity Linear 0-100% O₂

Zero Offset | <20µV

Temperature Range | -20°C to +50°C

Temperature Compensation | <2% variation from 0°C to 40°C

(see graph)

Differential Pressure Range 0-500mbar Max

Absolute Pressure Range | 500-2000mbar

Relative Humidity Range | 0 to 99% non-condensing | Cong Term Output Drift | <10% signal loss/year

Recommended Load Min 10KΩ

Resistor

Warranty Period | 12 month from date of despatch

N.B. All performance data is based on conditions at 20°C, 50%RH, and 1013mBar

NOTE

Molex header used in sensor is MOLEX 22-29-2031 Suggested mating parts are:

Molex 22-01-2035: 3-way housing

Molex 08-56-0110: crimp terminals

AO2 CiTiceL to be assembled into application 'finger tight' only

All tolerances ± 0.15mm unless otherwise stated

^{*} T_{99.5} response is equivalent to a change in concentration from 20.9% O₂ to 0.1% O₂

Oxygen CiTiceL® Specification



AO3 CiTiceL®

with jack socket connector

Performance Characteristics

Output | 9 - 13mV in Air

Range 0-100% O₂

Resolution 0.01% O₂

Expected Operating Life | 360000%O₂hrs at 20°C

286000%O₂hrs at 40°C

or 2 years in air at STP

T₉₀ Response Time <5 seconds T_{90 5} Response Time* <40 seconds

Signal in 100%O₂ | 100±1%

Linearity Linear 0-100% O₂

Zero Offset | <20µV

Temperature Range | -20°C to +50°C

Temperature Compensation | <2% variation from 0°C to 40°C

(see graph)

Differential Pressure Range 0-500mbar Max

Absolute Pressure Range | 500-2000mbar

Relative Humidity Range 0 to 99% non-condensing


Long Term Output Drift | <10% signal loss/year


Recommended Load | Min 10KΩ

Resistor

Warranty Period | 12 month from date of despatch

- * $T_{99.5}$ response is equivalent to a change in concentration from 20.9% O_2 to 0.1% O_2
- N.B. All performance data is based on conditions at 20°C, 50%RH, and 1013mBar





Doc. Ref.: a03.pmd Issue 1.7 Page 1 of 2 22nd October 2001